

Date/Time: 5/22/2020 8:15:41 AM
 Username: FSVC

UDP Report

Instrument #: CT1557(1)
 Version #: V4.72 A

Plasm-350002AMAX

Test Name	Version	Extraction Type	Master Mix Format					Ct Calculation	Lysis Heat Time
Notes	Channel Settings	475/520	530/565	585/630	630/665	680/715	Lysis Temperature	Sample Tip Height	
							Sample Volume	Wash Volume	
							Neutralization Volume	DNase Heat Time	
							Liquid Level Sensing Enabled		
Plasm-350002AMAX	3	ExK DNA-1 (Plasma/Serum) - 4-Snap	Type 4: Lyophilized MM with Rehydration Buffer					Call Ct at Threshold Crossing	20
--		Alias	Pfalc	Pviv	RNaseP	Pspg	37	1600	
		Gain	60	60	60	60	500	500	
		Threshold	200	200	200	200	--	--	
							Disabled		

Step Name: New Step
Profile Type: Hold **Cycles:** 1

Type	Time (s)	Temp (°C)	Detect
--	120	99	No

Step Name: New Step
Profile Type: 2 - Temperature **Cycles:** 45

Type	Time (s)	Temp (°C)	Detect
--	12	99	No
--	36.7	62	Yes



Performed By: _____ Date: _____
 Reviewed By: _____ Date: _____

Printed: 2020-05-22 08:15:41

Page: 1/1