

Date/Time: 5/22/2020 8:38:27 AM
 Username: FSVC

UDP Report

Instrument #: CT1557(1)
 Version #: V4.72 A

CMV 350053MAX

Test Name	Version	Extraction Type	Master Mix Format					Ct Calculation	Lysis Heat Time
Notes	Channel Settings	475/520	530/565	585/630	630/665	680/715	Lysis Temperature	Sample Tip Height	
							Sample Volume	Wash Volume	
							Neutralization Volume	DNase Heat Time	
							Liquid Level Sensing Enabled		
							CMV 350053MAX	2	ExK DNA-1 (Urine) - 4-Snap
--		Alias	--	--	CMV	--	SPC	62	1600
		Gain	--	--	60	--	60	937.5	500
		Threshold	--	--	200	--	200	--	--
								Disabled	

Step Name: New Step
Profile Type: Hold **Cycles:** 1

Type	Time (s)	Temp (°C)	Detect
--	300	99	No

Step Name: New Step
Profile Type: 2 - Temperature **Cycles:** 40

Type	Time (s)	Temp (°C)	Detect
--	14.2	99	No
--	44.6	62	Yes



Performed By: _____ Date: _____
 Reviewed By: _____ Date: _____

Printed: 2020-05-22 08:38:27