

Date/Time: 5/22/2020 8:13:43 AM
 Username: FSVC

UDP Report

Instrument #: CT1557(1)
 Version #: V4.72 A

PanAsp350058AMAX

Test Name	Version	Extraction Type	Master Mix Format				Ct Calculation		Lysis Heat Time
			475/520	530/565	585/630	630/665	680/715	Lysis Temperature	Sample Tip Height
Notes	Channel Settings						Sample Volume	Wash Volume	
							Neutralization Volume	DNase Heat Time	
							Liquid Level Sensing Enabled		
PanAsp350058AMAX	2	ExK DNA-1 (Urine) - 4-Snap	Type 4: Lyophilized MM with Rehydration Buffer				Call Ct at Threshold Crossing		9
--		Alias	--	--	Pan Asp	--	SPC	62	1600
		Gain	--	--	60	--	60	500	500
		Threshold	--	--	400	--	200	--	--
		Ct Min	--	--	--	--	--	--	--
		Ct Max	--	--	32	--	--	--	--
								Disabled	

Step Name: New Step			
Profile Type: Hold		Cycles: 1	
Type	Time (s)	Temp (°C)	Detect
--	300	99	No

Step Name: New Step			
Profile Type: 2 - Temperature		Cycles: 40	
Type	Time (s)	Temp (°C)	Detect
--	14	99	No
--	46.6	64	Yes



Performed By: _____ Date: _____
 Reviewed By: _____ Date: _____

Printed: 2020-05-22 08:13:44

Page: 1/1