

Date/Time: 1/30/2020 11:34:56 PM  
 Username: ADMIN

# UDP Report

Instrument #: MX0000(1)  
 Version #: V4.72 A

## 350068CMAX

Test Name	Version	Extraction Type	Master Mix Format					Ct Calculation		Lysis Heat Time
			475/520	530/565	585/630	630/665	680/715	Lysis Temperature	Sample Tip Height	
Notes	Channel Settings							Sample Volume	Wash Volume	
								Neutralization Volume	DNase Heat Time	
								Liquid Level Sensing Enabled		
350068CMAX	1	ExK DNA-1 (Urine) - 4-Snap	Type 4: Lyophilized MM with Rehydration Buffer					Call Ct at Threshold Crossing		20
		Alias	Cpsit	Cpne	Legspp	Mpne	SPC	62	1600	
		Gain	60	60	60	60	60	700	500	
		Threshold	200	200	200	200	200	--	--	
--		Excitation Channel 475/520	--	1	0	0	0			
		Excitation Channel 530/565	2.5	--	1	0	0			
		Excitation Channel 585/630	0	0	--	3	0			
		Excitation Channel 630/665	0	0	4	--	0			
		Excitation Channel 680/715	0	0	0	0	--			

Disabled

**Step Name:** New Step  
**Profile Type:** Hold **Cycles:** 1

Type	Time (s)	Temp (°C)	Detect
--	300	99	No

**Step Name:** New Step  
**Profile Type:** 3 - Temperature **Cycles:** 40

Type	Time (s)	Temp (°C)	Detect
--	10	99	No
--	25	62	Yes
--	13.7	72	No



Performed By: \_\_\_\_\_ Date: \_\_\_\_\_  
 Reviewed By: \_\_\_\_\_ Date: \_\_\_\_\_

Printed: 2020-01-30 23:34:57

Page: 1/1