

UDP Report

400003CMAX S4

Test Name	Version	Extraction Type	Master Mix Format				Ct Calculation	Lysis Heat Time
Notes	Channel Settings	475/520	530/565	585/630	630/665	680/715	Lysis Temperature	Sample Tip Height
							Sample Volume	Wash Volume
							Neutralization Volume	DNase Heat Time
							Liquid Level Sensing Enabled	
400003CMAX S4	1	ExK DNA-1 (Urine) - 4-Snap	Type 4: Lyophilized MM with Rehydration Buffer				Call Ct at Threshold Crossing	9
		Alias	Mgen	Urea	Mhom	Uparv	SPC	62
		Gain	60	60	60	60	60	500
		Threshold	200	200	200	200	200	--
		Excitation Channel 475/520	--	0	0	0	0	--
--		Excitation Channel 530/565	0	--	0	0	0	--
		Excitation Channel 585/630	0	0	--	5	0	--
		Excitation Channel 630/665	0	0	0	--	0	--
		Excitation Channel 680/715	0	0	0	0	--	--

Enabled

Step Name: New Step
Profile Type: Hold **Cycles:** 1

Type	Time (s)	Temp (°C)	Detect
--	300	99	No

Step Name: New Step
Profile Type: 2 - Temperature **Cycles:** 40

Type	Time (s)	Temp (°C)	Detect
--	14.4	99	No
--	44.4	62	Yes

	Mgen			
Mgen - PCR	✓	✓	✗	✗
SPC - PCR	✓	✗	✓	✗
Results	➕ POS	➕ POS	➖ NEG	⚠ UNR

	Urea			
Urea - PCR	✓	✓	✗	✗
SPC - PCR	✓	✗	✓	✗
Results	➕ POS	➕ POS	➖ NEG	⚠ UNR

	Mhom			
Mhom - PCR	✓	✓	✗	✗
SPC - PCR	✓	✗	✓	✗
Results	➕ POS	➕ POS	➖ NEG	⚠ UNR

	Uparv			
Uparv - PCR	✓	✓	✗	✗
SPC - PCR	✓	✗	✓	✗
Results	➕ POS	➕ POS	➖ NEG	⚠ UNR



Performed By: _____

Date: _____

Printed: 2018-10-18 13:47:49

Reviewed By: _____

Date: _____

Page: 1/1