

Date/Time: **10/18/2018 5:10:19 PM**
 Username: **ADMIN**

UDP Report

Instrument #: **CR0009(1)**
 Version #: **V4.70 A**

400006CMAX S4

Test Name	Version	Extraction Type	Master Mix Format				Ct Calculation	Lysis Heat Time
Notes	Channel Settings	475/520	530/565	585/630	630/665	680/715	Lysis Temperature	Sample Tip Height
							Sample Volume	Wash Volume
							Neutralization Volume	DNase Heat Time
							Liquid Level Sensing Enabled	
400006CMAX S4	1 ExK DNA-2 - 4-Snap	Type 4: Lyophilized MM with Rehydration Buffer				Call Ct at Threshold Crossing	9	
--	Alias	GBS	List	Ecoli	--	SPC	62	1600
	Gain	60	60	40	--	60	937.5	500
	Threshold	200	250	400	--	200	--	--
	Excitation Channel 475/520	--	3	0	0	0		
	Excitation Channel 530/565	3	--	0	0	0		
	Excitation Channel 585/630	0	0	--	0	0		
	Excitation Channel 630/665	0	0	0	--	0		
	Excitation Channel 680/715	0	0	0	0	--		

Disabled

Step Name: New Step
Profile Type: Hold **Cycles:** 1

Type	Time (s)	Temp (°C)	Detect
--	300	99	No

Step Name: New Step
Profile Type: 3 - Temperature **Cycles:** 40

Type	Time (s)	Temp (°C)	Detect
--	12	99	No
--	24.6	60	Yes
--	10	72	No

	GBS			
GBS - PCR	✓	✓	✗	✗
SPC - PCR	✓	✗	✓	✗
Results	➕ POS	➕ POS	➖ NEG	⚠ UNR

	List			
List - PCR	✓	✓	✗	✗
SPC - PCR	✓	✗	✓	✗
Results	➕ POS	➕ POS	➖ NEG	⚠ UNR

	Ecoli			
Ecoli - PCR	✓	✓	✗	✗
SPC - PCR	✓	✗	✓	✗
Results	➕ POS	➕ POS	➖ NEG	⚠ UNR



Performed By: _____ Date: _____
 Reviewed By: _____ Date: _____

Printed: 2018-10-18 17:10:19