

Date/Time: **11/27/2018 7:54:41 AM**  
 Username: **ADMIN**

# UDP Report

Instrument #: **CT0156(1)**  
 Version #: **V4.70 A**

## 400010WMAX S4

Test Name	Version	Extraction Type	Master Mix Format				Ct Calculation	Lysis Heat Time	
Notes	Channel Settings	475/520	530/565	585/630	630/665	680/715	Lysis Temperature	Sample Tip Height	
							Sample Volume	Wash Volume	
							Neutralization Volume	DNase Heat Time	
							Liquid Level Sensing Enabled		
400010WMAX S4	1	ExK DNA-1 (Urine) - 4-Snap	Type 4: Lyophilized MM with Rehydration Buffer				Call Ct at Threshold Crossing	20	
--		Alias	MTB	--	--	--	SPC	62	1600
--		Gain	60	--	--	--	60	700	500
--		Threshold	100	--	--	--	100	--	--
							Enabled		

**Step Name:** New Step  
**Profile Type:** Hold **Cycles:** 1

Type	Time (s)	Temp (°C)	Detect
--	120	99	No

**Step Name:** New Step  
**Profile Type:** 2 - Temperature **Cycles:** 45

Type	Time (s)	Temp (°C)	Detect
--	12	99	No
--	46.7	62	Yes

MTB				
MTB - PCR	✓	✓	✗	✗
SPC - PCR	✓	✗	✓	✗
Results	⊕ POS	⊕ POS	⊖ NEG	⚠ UNR



Performed By: \_\_\_\_\_ Date: \_\_\_\_\_  
 Reviewed By: \_\_\_\_\_ Date: \_\_\_\_\_

Printed: 2018-11-27 07:54:42