

Date/Time: 10/19/2018 4:35:40 PM
 Username: ADMIN

UDP Report

Instrument #: CT0156(1)
 Version #: V4.70 A

400012CMAX S4

Test Name	Version	Extraction Type	Master Mix Format				Ct Calculation	Lysis Heat Time
Notes	Channel Settings	475/520	530/565	585/630	630/665	680/715	Lysis Temperature	Sample Tip Height
							Sample Volume	Wash Volume
							Neutralization Volume	DNase Heat Time
							Liquid Level Sensing Enabled	
							400012CMAX S4	1 ExK DNA-3 - 4-Snap
	Alias	Glab	Alb	TropPar	Krus	SPC	62	
	Gain	80	80	80	80	80	700	
	Threshold	200	200	300	200	200	--	
	Excitation Channel 475/520	--	0	0	0	0	--	
--	Excitation Channel 530/565	5	--	5	0	0		
	Excitation Channel 585/630	0	0	--	5	0		
	Excitation Channel 630/665	0	0	5	--	0		
	Excitation Channel 680/715	0	0	0	0	--		

Enabled

Step Name: New Step
Profile Type: Hold **Cycles:** 1

Type	Time (s)	Temp (°C)	Detect
--	300	99	No

Step Name: New Step
Profile Type: 2 - Temperature **Cycles:** 40

Type	Time (s)	Temp (°C)	Detect
--	15	99	No
--	47.2	66	Yes

	Glab			
Glab - PCR	✓	✓	✗	✗
SPC - PCR	✓	✗	✓	✗
Results	⊕ POS	⊕ POS	⊖ NEG	⚠ UNR

	Alb			
Alb - PCR	✓	✓	✗	✗
SPC - PCR	✓	✗	✓	✗
Results	⊕ POS	⊕ POS	⊖ NEG	⚠ UNR

	TropPar			
TropPar - PCR	✓	✓	✗	✗
SPC - PCR	✓	✗	✓	✗
Results	⊕ POS	⊕ POS	⊖ NEG	⚠ UNR

	krus			
Krus - PCR	✓	✓	✗	✗
SPC - PCR	✓	✗	✓	✗
Results	⊕ POS	⊕ POS	⊖ NEG	⚠ UNR



Performed By: _____

Date: _____

Printed: 2018-10-19 16:35:41

Reviewed By: _____

Date: _____

Page: 1/1