

Date/Time: **10/23/2018**  
**10:07:45 AM**  
 Username: **ADMIN**

# UDP Report

Instrument #: **CR0009(1)**  
 Version #: **V4.70 A**

## 400037CMAX S4

Test Name	Version	Extraction Type	Master Mix Format				Ct Calculation		Lysis Heat Time
			475/520	530/565	585/630	630/665	680/715	Lysis Temperature	Sample Tip Height
Notes	Channel Settings						Sample Volume	Wash Volume	
							Neutralization Volume	DNase Heat Time	
							Liquid Level Sensing Enabled		
400037CMAX S4	1	ExK DNA-2 - 4-Snap	Type 4: Lyophilized MM with Rehydration Buffer				Call Ct at Threshold Crossing		9
--		Alias	vanA	--	vanB	--	SPC	62	1600
		Gain	60	--	60	--	60	937.5	500
		Threshold	200	--	200	--	200	--	--
								Enabled	

**Step Name:** New Step  
**Profile Type:** Hold **Cycles:** 1

Type	Time (s)	Temp (°C)	Detect
--	300	99	No

**Step Name:** New Step  
**Profile Type:** 2 - Temperature **Cycles:** 40

Type	Time (s)	Temp (°C)	Detect
--	12.8	99	No
--	46	62	Yes

VanA				
vanA - PCR	✓	✓	✗	✗
SPC - PCR	✓	✗	✓	✗
Results	⊕ POS	⊕ POS	⊖ NEG	⚠ UNR

VanB				
vanB - PCR	✓	✓	✗	✗
SPC - PCR	✓	✗	✓	✗
Results	⊕ POS	⊕ POS	⊖ NEG	⚠ UNR



Performed By: \_\_\_\_\_

Date: \_\_\_\_\_

Printed: 2018-10-23 10:07:46

Reviewed By: \_\_\_\_\_

Date: \_\_\_\_\_

Page: 1/1