

Xfree™ COVID-19 Direct RT-PCR

Workflows for Direct and Extracted Samples

Quick guide

High-Throughput Workflow

80 minutes for 96 to 384 samples

Validated platforms: Applied Biosystems™ 7500 Fast Dx and QuantStudio™ 5
Bio-Rad CFX96 Touch™ and CFX384 Touch™

Kit Component

Xfree™ Sample-Ready™ Tube



Supplied by BioGX


Molecular Grade Water



Purchased separately


Materials Required

96/384-well PCR Plate or PCR Strip Tubes



Purchased separately

Optical Plate Seals or Tube Caps



Purchased separately


Direct Sample Workflow

Validated Transport Media:

Saline
UTM®
UVT
VTM
ESwab™



Rehydrate with 400 μ L of molecular grade water. Dispense Xfree™ reagent (24 x 15 μ L) into multi-well plate.



Add 5 μ L of patient sample, pipette up and down once and apply the optical seal or tube caps.



Pulse spin, load into real-time PCR instrument & start run protocol specific for direct sample.

Extracted Sample Workflow

Validated Transport Media:

Saline
UTM®
UVT
VTM
ESwab™



Rehydrate with 400 μ L of molecular grade water. Dispense Xfree™ reagent (36 x 10 μ L) into multi-well plate.



Add 5 μ L of extracted RNA, pipette up and down once and apply the optical seal or tube caps.



Pulse spin, load into real-time PCR instrument & start run protocol specific for extracted sample.

